

## DATASHEET


### Lenti-mKate2-Neo virus

Cat. VSL-0033


Titer:	~10 <sup>6</sup> IU/ml
Medium:	DMEM, 10% FBS
Volume:	10 ml (10x1ml)
Selection marker:	Neomycin (G418)
Storage:	-70°C
How to use:	Thaw the recombinant lentivirus supernatant in a 37°C water bath; remove it from the bath immediately when thawed.

**Description:** Ready-to-use virus expressing far-red fluorescent protein mKate2 can be used as a control in experiments with lentivirus transductions of target cells or for the generation of stable cell lines.

mKate2 is the next generation of monomeric far-red fluorescent protein TagFP635 (mKate) [Shcherbo et al., 2007; Shcherbo et al., 2009]. Lenti-Mkate2 virus is VSV-G pseudotyped lentivirus that can be used as a control for infection of target cells. It was generated by co-transfection of 293T packaging cells with pLVmKate2 plasmid (Fig.1) and Packing Mix. Virus-containing supernatant has been harvested 48-72 hours post-transfection, centrifuged at 3000rpm for 15 minutes at +4°C to pellet debris, and then filtered through 0.45µm PVDF syringe.



**Fig.1. Map of pLVmKate2-Neo vector**



**Fig.2. 293 cells infected with Lenti-mKate2 virus**

Notice to Purchaser: All products containing mKate2 are intended for research use only. The Products are covered by U.S. Pat. # 7,417,131 and other Evrogen Patents and/or Patent applications pending. By use of these Products, you accept the terms and conditions of the applicable Limited Use Label License.

***For research use only***

9700 Great Seneca Hwy, Suite 247  
Rockville, MD 20850  
Tel: 800-475-2812; Fax: 240-363-0739